

Chapter 1

WLM class

WLM is used to control a High Finesse WS6 or WS7 wavemeter. Only one WLM object can be instantiated. The company provides a C-library of functions to configure and control the wavemeter, hence the WLM object loads this library ('C:\Windows\System32\wlmData.dll') and the corresponding C-header ('wlmData.hml') in the constructor function. C-functions from the library are called within the matlab functions described below. A complete summary of the C functionality can be found in the wavemeter manual, provided by High Finesse.

MATLAB EXAMPLE

```
handle = WLM.getInstance()  
output = handle.[anyFunction];
```

Properties

- *WL*
Last measured wavelength in nanometers.
- *freq*
Last measured frequency in THz.
- *pool*
Handle to parallel worker pool object. Used for toggling the wavemeter

channel in parallel with other matlab scripts for simultaneous wavemeter locking of multiple lasers.

- *toggleing*
String that shows whether toggleing of the wavemeter channels is 'on' or 'off'.
- *parf*
Handle to parallel function, executed by the parallel worker pool *pool*.
- *ar* (private)
Struct containing the argument name strings from WLM header.
- *active_channel* (private)
Channel that is currently active.
- *num_channels* (private)
Total number of switchbox channels.

Functions

handle = **getInstance()**

Instantiates a WLM object and puts it under the name *handle* (or any other name that is put before the = sign).

obj = **WLM()** (private)

Constructor that is called by `getInstance()`. Here, the C-header is called, its parameters are stored in private property *ar*. Also, the C-library with functions to communicate with the wavemeter is loaded. To enable parallel execution of functions, the library is also loaded into the workers in *pool*.

[*WL* , *freq*] = **ReadWLM(what , channel)** (private)

Switches the wavemeter to switchbox *channel* and reads frequency and/or wavelength, depending on whether *what* is 'WLM', 'freq' or 'both' (cap-sensitive!). Called by `GetWL()`, `GetFreq()` and `GetBoth()`.

$pool = \mathbf{CreateParPool}()$ (private)

Checks if there is a parallel worker pool available and stores it in output $pool$. If not available, it creates a new one. This might take a minute.

$\mathbf{ParToggle}(channels , time)$ (private)

Starts an infinite loop that switches to the next channel in vector $channels$ every $time$ seconds. Called as a parallel function by `ToggleChannels()`.

$\mathbf{ToggleChannels}(channels , time , onOffStr)$

Is called by the user to start toggling between channels for simultaneous locking of multiple lasers. If string $onOffStr$ is 'on', function `ParToggle(channels , time)` is started as a job in the parallel worker pool. A handle to this job is stored in property $parf$. If $onOffStr$ is 'off', the job stored in $parf$ is canceled and the toggling stops.

$\mathbf{SwitchToChannel}(channel)$

Changes the active switchbox channel to $channel$, then waits for 0.5s.

$freq = \mathbf{GetFreq}(channel)$

Calls `ReadWLM('freq' , channel)` to measure the frequency for $channel$.

$WL = \mathbf{GetWL}(channel)$

Calls `ReadWLM('WLM' , channel)` to measure the wavelength for $channel$.

$[WL , freq] = \mathbf{GetBoth}(channel)$

Calls `ReadWLM('both' , channel)` to measure both frequency and wavelength for $channel$.

channel = **GetChannel()**

Returns the property *active_channel*.